

Metabolite profiling of blood specimens using a gas chromatography time of flight mass spectrometry (GCTOFMS)

West Coast Metabolomics Center, University of California, Davis, USA
www.metabolomics.ucdavis.edu

Summary

This SOP describes the instructions for generating metabolomics datasets for blood specimens using a GCTOFMS instrument. It includes the extraction of metabolites, the acquisition of mass spectrometry data and the data processing and normalization to generate final data matrices.

Method

Extraction of metabolites from blood Plasma specimens

1. References:

- Fiehn O, Kind T (2006) Metabolite profiling in blood plasma. In: Metabolomics: Methods and Protocols. Weckwerth W (ed.), Humana Press, Totowa NJ (in press)

2. Starting material:

- Blood plasma/serum: 30 µl sample volume (validated for EDTA and heparin plasma, and serum)

3. Equipment:

- Centrifuge (Eppendorf 5415 D)
- Calibrated pipettes 1-200µl and 100-1000µl
- Eppendorf tubes 2ml, uncoloured (Cat.No. 022363204)
- ThermoElectron Neslab RTE 740 cooling bath at -20°C
- MiniVortexer (VWR)
- Orbital Mixing Chilling/Heating Plate (Torrey Pines Scientific Instruments)
- Speed vacuum concentration system (Labconco Centrивap cold trap)

4. Chemicals


- Acetonitrile, LCMS grade (JT Baker; Cat. No.9829-02)
- Isopropanol, HPLC (JT Baker; Cat. No. 9095-02)
- Crushed ice
- pH paper 5-10 (EMD Chem. Inc.)
- Nitrogen line with pipette tip
- 18 MΩ pure water (Millipore)

5. Preparation of extraction mix and material before experiment:

1. Switch on bath to pre-cool at -20°C (±2°C validity temperature range)
2. Check pH of acetonitrile and isopropanol (pH7) using wetted pH paper
3. Make the extraction solution by mixing acetonitrile, isopropanol and water in proportions 3 : 3 : 2
4. Rinse the extraction solution mix for 5 min with nitrogen. Make sure that the nitrogen line was flushed out of air before using it for degassing the extraction solvent solution

6. Sample preparation:

1. Vortex the plasma/serum samples for 10s to obtain a homogenized sample using the MiniVortexer.
2. Aliquot 15-30ul and add 1mL extraction solution. The extraction solution has to be pre-chilled using the ThermoElectron Neslab RTE 740 cooling bath set to -20°C.



3. Vortex the sample for 10s and shake for 5min at 4°C using the Orbital Mixing Chilling/Heating plate. If you are using more than one sample, keep the rest of the samples on ice (chilled at <0°C with NaCl).
4. Centrifuge samples for 2min at 14000 rcf using the centrifuge Eppendorf 5415 D.
5. Aliquot two 500 µL portions of the supernatant. One for analysis and one for backup. Store one aliquot in the -20°C freezer as a backup.
6. Evaporate one 500 µL aliquot of the sample in the Labconco Centrивap cold trap concentrator to complete dryness.
7. The dried aliquot is then re-suspended with 500 µL 50% acetonitrile (degassed as given above).
8. Vortex the samples for 10s using the MiniVortexer VWR.
9. Centrifuge for 2min at 14000 rcf using the centrifuge Eppendorf 5415 D.
10. Remove supernatant to a new Eppendorf tube.
11. Evaporate the supernatant to dryness in the Labconco Centrивap cold trap concentrator.
12. Submit to derivatization.

6. Problems

To prevent contamination disposable material is used. Control pH from extraction mix.

7. Quality assurance

For each sequence of sample extractions, perform one blank negative control extraction by applying the total procedure (i.e. all materials and plastic ware) without biological sample.

8. Disposal of waste

Collect all chemicals in appropriate bottles and follow the disposal rules.

Acquisition of GCMS raw data

Instruments:

**Gerstel CIS4 –with dual MPS Injector/
Agilent 6890 GC- Leco Pegasus III TOF MS**

Injector conditions:

Agilent 6890 GC is equipped with a Gerstel automatic liner exchange system (ALEX) that includes a multipurpose sample (MPS2) dual rail, and a Gerstel CIS cold injection system (Gerstel, Muehlheim, Germany) with temperature program as follows: 50°C to 275°C final temperature at a rate of 12 °C/s and hold for 3 minutes. Injection volume is 0.5 µl with 10 µl/s injection speed on a splitless injector with purge time of 25 seconds. Liner (Gerstel #011711-010-00) is changed after every 10 samples, (using the Maestro1 Gerstel software vs. 1.1.4.18). Before and after each injection, the 10 µl injection syringe is washed three times with 10 µl ethyl acetate.

Gas Chromatography conditions:

A 30 m long, 0.25 mm i.d. Rtx-5Sil MS column (0.25 µm 95% dimethyl 5% diphenyl polysiloxane film) with additional 10 m integrated guard column is used (Restek, Bellefonte PA). 99.9999% pure Helium with built-in purifier (Airgas, Radnor PA) is set at constant flow of 1 ml/min. The oven temperature is held constant at 50°C for 1 min and then ramped at 20°C/min to 330°C at which it is held constant for 5 min.

Mass spectrometer settings:

A Leco Pegasus IV time of flight mass spectrometer is controlled by the Leco ChromaTOF software vs. 2.32 (St. Joseph, MI). The transfer line temperature between gas chromatograph and mass spectrometer is set to 280°C. Electron impact ionization at 70V is employed with an ion source temperature of 250°C. Acquisition rate is 17 spectra/second, with a scan mass range of 85-500 Da.

Data processing and normalization :

Data processing Raw data files are preprocessed directly after data acquisition and stored as ChromaTOF-specific *.peg files, as generic *.txt result files and additionally as generic ANDI MS *.cdf files. ChromaTOF vs. 2.32 is used for data preprocessing without smoothing, 3 s peak width, baseline subtraction just above the noise level, and automatic mass spectral deconvolution and peak detection at signal/noise levels of 5:1 throughout the chromatogram. Apex masses are reported for use in the BinBase algorithm. Result *.txt files are exported to a data server with absolute spectra intensities and further processed by a filtering algorithm implemented in the metabolomics BinBase database. The BinBase algorithm (rtx5) used the settings: validity of chromatogram (10^7 counts s⁻¹), unbiased retention index marker detection (MS similarity > 800, validity of intensity range for high m/z marker ions), retention index calculation by 5th order polynomial regression. Spectra are cut to 5% base peak abundance and matched to database entries from most to least abundant spectra using the following matching filters: retention index window ±2,000 units (equivalent to about ±2 s retention time), validation of unique ions and apex masses (unique ion must be included in apexing masses and present at >3% of base peak abundance), mass spectrum similarity must fit criteria dependent on peak purity and signal/noise ratios and a final isomer filter. Failed spectra are automatically entered as new database entries if s/n > 25, purity 80%. All thresholds reflect settings for ChromaTOF v. 4.0. Quantification is reported as peak height using the unique ion as default, unless a different quantification ion is manually set in the BinBase administration software BinView. A quantification report table is produced for all database entries that are positively detected in more than 10% of the samples of a study design class (as defined in the miniX database) for unidentified metabolites. A subsequent post-processing module is employed to automatically replace missing values from the *.cdf files. Replaced values are labeled as 'low confidence' by color coding, and for each metabolite, the number of high-confidence peak detections is recorded as well as the ratio of the average height of replaced values to high-confidence peak detections. These ratios and numbers are used for manual curation of automatic report data sets to data sets released for submission.

Data reporting


Data are reported including metadata, see example below.

							Subject ID	223913	157819	124940
							Local code	A0050702A	A0125621A	142363
							Vial Barcode	1RAR7	1GZR9	1AN1N
							Date received	4-Dec-12	4-Dec-12	4-Dec-12
							Date of evaluation	3/17/2013	3/17/2013	3/17/2013
							Sample Status			
							REVISION			
							Comments	GCTOF MS_s	GCTOF MS_sa	GCTOF MS_s
							Acq. Time	5:18:33 AM	5:43:50 AM	6:09:05 AM
							Data File Name	130328cmss	130328cmss	130328cmss
							miniXid	118072	118073	118074
BinBase i	BinBase name	ret.index	quant m/z	mass spec	InChI key	KEGG id	PubChem id			
14441	z C30 FAME internal stand	1113100	87	82:386.0 83:	BIRUBGLRQLAEFF-UHFFF	n/a	12400	16026	15203	18096
14378	z C28 FAME internal stand	1061700	87	82:1635.0 83	ZKHOYAKAFALNQD-UHFF	n/a	41518	39317	228	11145
14367	z C26 FAME internal stand	1006900	87	82:1915.0 83	VHUJBYFFWDLNM-UHFF	n/a	22048	32809	30571	35507
14373	z C24 FAME internal stand	948820	87	82:1153.0 83	XUDJZDNUVZHSKZ-UHFFF	n/a	75546	43836	43163	48731
14350	z C22 FAME internal stand	886620	87	82:3004.0 83	QSQLTHHMFHEFY-UHFF	n/a	13584	53566	51781	58740
14338	z C20 FAME internal stand	819620	87	82:5074.0 83	QGBRLVONZXHAKJ-UHFF	n/a	14259	57778	58877	63250
14344	z C18 FAME internal stand	747420	87	82:1871.0 83	HPEUJPOZXNMSJ-UHFFFA	n/a	8201	53466	52542	57091
14328	z C16 FAME internal stand	668720	87	82:7874.0 83	FLIACVVOZYBSBS-UHFFF	C16995	8181	170150	166843	186189
14330	z C14 FAME internal stand	582620	87	82:6846.0 83	ZAZKJZBWRNLDLS-UHFF	n/a	31284	117029	113463	131850
15538	z C12 FAME internal stand	487220	87	87:89125.0 1	UQDUPQYQJQYHUI-UHFFF	n/a	8139	147184	142864	169095
14348	z C10 FAME internal stand	381020	87	82:2900.0 83	YRHYCMZPEVDGFQ-UHFF	n/a	8050	131518	127669	150456
14356	z C09 FAME internal stand	323120	87	82:1461.0 83	IJXHLVMUNBOGRU-UHFF	n/a	15606	119374	119013	137532
14391	z C08 FAME internal stand	262320	87	82:816.0 83:	JGHZJRVDZXSNNQ-UHFFF	n/a	8091	79355	87424	91461
231968	xylose	544673	103	85:5068.0 86	PYMYPHUHKUWMLA-VPE	C02205	644160	19097	1462	1642
368041	xylose	566570	217	86:588.0 87:	HEBKCHPVOIAQTA-NGQZ	C00379	6912	239	168	121
203224	xanthine	702391	353	85:1361.0 86	LRFTTYWOQMYALW-UHF	C00385	1188	64	26	60
199605	valine	313224	144	85:48.0 86:1	KZSNJWFQEVHDMF-BYPY	C00183	6287	108089	120432	133290
213127	uridine	856953	258	85:2472.0 86	DRTQHJPVMGBCUF-XVFC	C00299	6029	256	63	60
304993	uric acid	730534	441	85:1183.0 86	LEHOTFFKMJEONL-UHFF	C00366	1175	16528	13715	7856
224322	urea	337230	171	87:1643.0 89	XSQKJJFZCRTK-UHFFFA	C00086	1176	334000	281631	313888

The 'BinBase identifier column' denotes the unique identifier for the GCTOFMS platform. It is given for both identified and unidentified metabolites in the same manner. The 'BinBase name' denotes the name of the metabolite, if the peak has been identified. A chemical name is not a unique identifier. We use names recognized by biologists instead of IUPAC nomenclature. If a compound is identified, it has a name, and external database identifiers such as InChI key, PubChem ID and KEGG ID. If a compound is unknown, the name is the same as given in the 'identifier column'.

The 'retention index' column details the target retention index in the BinBase database system. The 'quant m/z' column details the m/z value that was used to quantify the peak height of a BinBase entry. The 'mass spec' column details the complete mass spectrum of the metabolite given as m/z: intensity values, separated by spaces. The 'InChI key' identifier gives the unique chemical identifier defined by the IUPAC and NIST consortia. The 'KEGG' identifier gives the unique identifier associated with an identified metabolite in the community database KEGG LIGAND DB. The 'PubChem' column denotes the unique identifier of a metabolite in the PubChem database. The 'internal standard' addition within the BinBase name clarifies if a specific chemical has been added into the extraction solvent as internal standard. These internal standards serve as retention time alignment markers, for quality control purposes and for quantification corrections. Row metadata that are requested by a specific consortium are labeled in blue. Consortium 'subject ID', 'local ID', 'vial barcode' detail information given by a specific consortium. The row 'date received' is the date when samples were received in the metabolomics laboratory. The row 'date of evaluation' is the data of data acquisition, as given by the machine logbook. The row 'sample status' uses the consortium's sample status code if samples have errors. The consortium sample status code does not give a code when data acquisition occurred without problems. If a consortium does not use an authorized error code dictionary, plain text is given for errors. The row 'revision' details if data processing yields a new data sheet. Data revisions may be needed when new algorithms have been tested, validated and deployed that might yield better raw data analyses than prior submissions. By default, therefore, data revisions replace the (less valid) prior data submissions. However, data revisions may also indicate a different form of data treatment, e.g. data normalizations (see below). In this case, the 'revision' would indicate the type of normalization. Any information in the row 'revision' will have a date stamp when the revision was conducted in the form of MMDDYY. The 'comments' row gives comments about the platform and type of

Rev Jan 31 2011



sample. A sample is given as “sample” in comparison to e.g. a quality control or a blank injection. The ‘Acq.Date-Time’ row details the acquisition time when the data acquisition was completed. The ‘Data File Name’ row denotes the name of the raw data file. Raw data files are secured at the NIH Metabolomics database, www.metabolomicsworkbench.org. Data file names are dictated by the laboratory’s information and management system when the sequence starts running. GCTOF raw file names from the Leco instrumentation end with .peg (this ending is not given in the file name, but is found in the database repositories). In case a sample will need to be reinjected, the file name will change from e.g. 130328cmssa40_1.peg to 130328cmssa40_2.d for the second injection, 130328cmssa40_3.d for the third injection and subsequent injections. The file name itself denotes YYMMDD then the ‘machine used for data acquisition’ (here: c; we have four GCTOF MS machines a-d), ‘person who operated the machine’ (here: ms for Mimi Swe), ‘sa’ for sample (instead of e.g. ‘qc’ for a quality control or ‘bl’ for a blank sample), followed by the sequence number (here: the 40th sample within the sample sequence). The ‘miniX’ row shows the unique sample identifier in the Fiehnlab miniX laboratory information management system.

The actual data are given as peak heights for the quantification ion (mz value) at the specific retention index. We give peak heights instead of peak areas because peak heights are more precise for low abundant metabolites than peak areas, due to the larger influence of baseline determinations on areas compared to peak heights. Also, overlapping (co-eluting) ions or peaks are harder to deconvolute in terms of precise determinations of peak areas than peak heights. Such data files are then called ‘raw results data’ in comparison to the raw data file produced during data acquisition (see ‘data file name’). The worksheets are called ‘Height’.

Raw results data need to be normalized to reduce the impact of between-series drifts of instrument sensitivity, caused by machine maintenance, aging and tuning parameters. Such normalization data sets are called ‘norm data’ worksheets.

There are many different types of normalizations in the scientific literature. We usually provide first a variant of a ‘vector normalization’ in which we calculate the sum of all peak heights for all identified metabolites (but not the unknowns!) for each sample. We call such peak-sums “mTIC” in analogy to the term TIC used in mass spectrometry (for ‘total ion chromatogram’), but with the notification “mTIC” to indicate that we only use genuine metabolites (identified compounds) in order to avoid using potential non-biological artifacts for the biological normalizations, such as column bleed, plasticizers or other contaminants.

Subsequently, we determine if the mTIC averages are significantly different between treatment groups or cohorts. If these averages indeed are different by $p < 0.05$, data will be normalized to the average mTIC of each group. If averages between treatment groups or cohorts are not different, or if treatment relations to groups are kept blinded, data will be normalized to the total average mTIC.

Following equation is then used for normalizations for metabolite i of sample j :

$$\text{metabolite}_{ij, \text{normalized}} = \frac{\text{metabolite}_{ij, \text{raw}}}{\text{mTIC}_j} \cdot \text{mTIC}_{\text{average}}$$

The worksheet is then called ‘norm mTIC’. Data are ‘relative semi-quantifications’, meaning they are normalized peak heights. Because the average mTIC will be different between series of analyses that are weeks or months apart (due to differences in machine sensitivity, tuning, maintenance status and other parameters), additional normalizations need to be performed. For this purpose, identical samples (‘QC samples’) must be analyzed multiple times in all series of data acquisitions. In fact, one must not exclude the possibility that even within a series of data acquisitions, a sensitivity shift or drift might occur. Hence, the following statistical analyses are suggested: (a) compute univariate statistics for mTIC values in batches within-series and between-series of data injections, using time/date stamps to find potential breaks during which machine downtime may have occurred. If there are no mTIC differences between such time/date stamp batches, calculate an overall mTIC covering all samples. (b) compute multivariate PCA plots for the , marking the potentially different samples of individual time/date stamp batches using different colors. If there is no apparent separation between PCA clusters of different colors, there is no large between-series effect and these PCA clusters can be treated as indistinguishable. If there is suspicion of hidden features that might be masked by overall variance analysis in PCA, supervised statistics by Partial Least Square regression models can unravel such between-series

differences. Once different clusters (i.e. series of undistinguishable QC samples) have been identified, correction factor models need to be developed that correct differences between those QC samples. Subsequently, these correction factors can be applied to the actual analytical samples to remove overt quantification differences that are not related to biological causes but solely due to analytical errors.

Such correction factor models can be computed in different ways, e.g. by unit-variance mean centering or by calculating simple offset vectors for each individual metabolite. The best way of such types of normalizations is being explored in the Fiehn laboratory. However, in any case, such correction models can only be developed if a sufficient number of QC samples have been included in the analytical sequences. For that reason, the Fiehn laboratory uses a suitable QC sample for every 11th injection. Such QC samples need to be as similar to the actual biological specimen as possible, e.g. generated by pool samples during extractions or by obtaining typical community standard samples (e.g. the NIST standard blood plasma, or commercial serum or plasma samples as needed).

If appropriate internal standards are used for absolute quantifications, the following equation could be used for peak height normalizations for metabolite i of sample j and internal standard k

$$\text{metabolite}_{ij, \text{normalized}} = \frac{\text{metabolite}_{ij, \text{raw}}}{\text{istd}_k \cdot \text{concentration istd}_k}$$

However, there are few universal or class-specific internal standards in GC-MS based analysis, because within each chemical class, metabolites may have drastically different calibration curves (sensitivity or 'response') based on a combination of injection, volatilization and stability and ionization response properties. As surrogate, external calibration standards could be used for specific (important) metabolites which, however, cannot be applied for unidentified compounds and which of course would not account for recovery during extraction procedures.

Filtering, normalization and preparation of GCTOF metabolomics dataset for the ADNI-I study

The analysis included 833 samples, 83QC and 4 NIST plasma. 143 known and 241 unknown were measured using gas chromatography and time of light mass spectrometry and binbase database processing. Data were acquired on two GCTOF instruments in a period of a month (7/6/2015 to 8/4/2015) in multiple batches. This report explains the different steps to diagnosis and fix the drifts in signals caused by batches on two machines and how we have generated a statistics-ready data matrix.

Step 1.

Detection of sample outliers: One Outlier detected from PCA score plot.




Figure 1. Because this outlier sample is not a QC, it was not removed it when normalizing.

Step 2. How to define the analysis batches:




Figure 2. Definition of analytical batches. There are two machines A and C. We define each machine as one batch. Furthermore, from the time plot above, we can see that there are five time intervals. For the study, batches were defined by longer breaks in analysis sequences on two machines, generating 7 batches in total. Thus we define batch as **machine * time interval**. As a result, there are in total seven batches, highlighted by different color on the following plot. Not all the batches of same size.




Figure 3. A and C are GCTOF machines and a-d are time intervals.

Step 3. Visual inspection of batch effect:

Since, FAME markers were added externally in the samples, they are the most suitable compounds to check for differences caused by machines and analysis time. Following is an example of batch effect on C30.FAME.internal.standard.




Figure 4. Effect of different batches on the signal of C30 fame marker. X-axis is peak intensity.

Step 4. LOESS (locally weighted scatterplot smoothing) and batch-median based signal correction

LOESS correction within a batch was used first then the median correction between bathes was utilized.

Median correction formula is following, for each sample,

After = before * ratio,

Where ratio is global median of QC divided by median of the batch median corresponding to that sample. As a result, QCs after normalization has the same median with QCs before normalization.




Figure 5. Upper panel: Raw signal intensities. Lower panel: Signal intensities after normalization.

Step 5. Optimization of LOESS parameters.

The method could fix the signal drift for a majority of compounds, but, for some compounds, the outlier of QC completely destroys the loess correction within bathes. LOESS regression is sensitive to outliers, and may introduce noise in data. For example the “A a” batch in the following plot, outliers make the loess line off the trend.




Figure 6. Upper panel : Effect of outliers in the LOESS regression line. Lower panel : after excluding the outlier in the LOESS regression.

Therefore before normalizing each batch, boxplots were used to find and exclude outliers (beyond 1.5 IQR) before loess correction in order to avoid that LOESS modeling is affected by outliers.

Step 6. Manual inspection and filtering of compounds:

After all compounds were normalized, scatter plot of each compound were examined, comparing patterns before and after normalizations. We found that some compounds had unexpected QC behavior (following is an example), with different trends of intensity patterns in the QC and the real samples, specifically for some compounds with very low intensities (<500). We made decisions on each compounds either delete those or keep un-normalized data as is. The criteria were compound

RSD, peak level and detection frequency.




Figure 7.

For five compounds, creatinine, maltose, maleimide, inosine and cysteine, un-normalized data was used as LOESS did not performed well for them. Also, for proline, normalization with total signals for known metabolites was performed. 135 known and 86 unknown metabolites passed the filtering. Rest of compounds were excluded from the data matrix. Targeted peak picking in mzmine was performed to re-create the pick-list.

Dataset Information

This methods document applies to the following dataset(s) available from the ADNI repository:

Dataset Name	Date Submitted
mx 232603 Rima Kaddurah-Daouk_ADNI_human plasma 06-2016_submit	06/01/2016

About the Authors

For more information please contact Dr Dinesh Kumar Barupal or Dr Oliver Fiehn by email at metabolomics@ucdavis.edu.

Notice: This document is presented by the author(s) as a service to ADNI data users. However, users should be aware that no formal review process has vetted this document and that ADNI cannot guarantee the accuracy or utility of this document.