

Longitudinal Analysis of Cerebrospinal Fluid Visinin-like protein-1

Courtney Sutphen, BS, Elizabeth Macy, BA, Jack Ladenson, PhD, Anne Fagan, PhD

| Contents | |
|----------|--|
| Page 1 | Summary and Method |
| Page 2 | Method (Cont.) and Dataset Information |
| Page 3 | References and About the Authors |


Summary

Recent preclinical biomarker studies in AD have shown that development of the pathological hallmarks of AD begins as early as 10-15 years before the development of cognitive decline¹⁻⁵. Clinicians, researchers, and pharmaceutical companies have begun to appreciate possible intervention in the preclinical stages of AD as a patient's best option for the prevention of dementia and, in the shorter term, as the best opportunity to test proposed disease-modifying drugs currently in development. Cerebrospinal fluid (CSF) Visinin-like protein-1 (VILIP-1) shows utility in distinguishing AD from cognitively normal (CN) controls and in some cases from other dementias^{6,7}. Importantly, it appears useful as a prognostic marker, particularly when used in combination with A β 42. Further work, including in-depth cross-sectional studies across multiple cohorts and especially longitudinal follow-up of CN subjects, is needed to improve the ability of biomarkers to contribute to AD research and clinical trials. Adding VILIP-1 and other novel biomarkers to a panel including A β 42, tau, and ptau181 may improve our diagnostic and prognostic capabilities. Better knowing the order of changes in biomarkers as preclinical AD progresses to symptomatic AD will enable clinical trials to select those participants who are most likely to benefit from drug intervention, as well as to reduce the number of participants needed to attain statistical relevance.

Method

Mouse anti-human VILIP-1, clone 3A8.1 and sheep anti-human VILIP-1 were provided to Singulex, Inc. (Alameda, CA) for development and manufacture of an Erenna® immunoassay assay kit (#03-0008-02) according to an agreement with Washington University. The proprietary custom order kit reagents included preparation of the monoclonal capture antibody for binding to Invitrogen (Carlsbad, CA) MyOne magnetic micro particles and Invitrogen Alexa fluor dye labeling of the sheep polyclonal detection antibody and all assay buffers.

A sandwich ELISA was developed using the Erenna® immunoassay system to measure VILIP-1 in CSF. Prior to the assay all samples were centrifuged (11,000 g x 3 minutes) to remove particulates. All assay steps were performed at room temperature unless otherwise indicated. A calibration curve was prepared using dilutions of recombinant human VILIP-1 ranging from 3.9 to 3000 pg/mL and , with each concentration assayed in triplicate. 15 μ L standards or CSF were combined with 135 μ L assay buffer and 50 μ L antibody coated micro



particles for measurement of VILIP-1 in CSF. The assay plate was incubated for two hours on a plate shaker set to 525 revolutions per minute. Micro particles were then magnetically separated and washed one time using an Agilent (Santa Clara, CA) Bravo Automated Liquid Handling Platform using Singulex wash buffer. Fluorescent dye labeled detection antibody prepared by Singulex (20 μ L per well) was added and incubated for one hour. After washing the magnetic micro particles 5 times, 20 μ L per well of Singulex elution buffer was added for 10 minutes to separate detection antibody from the micro particles. Eluted antibodies were then transferred with the Bravo instrument to a clean 384 well plate for reading in the Erenna® immunoassay system.

All samples (first freeze/thaw cycle) were run in duplicate on a single lot number, and within-person longitudinal samples were run on the same assay plate to reduce inter- and intra-plate variability. Quality control (QC) included analysis of three internal standard CSF pools run on each plate. Each standard was averaged over all assays; plates with 2 or more standards greater than two standard deviations (SD) from the mean were reanalyzed. Further, any individual sample with a Coefficient of Variation (% CV) greater than 25% was reanalyzed. Samples that failed quality control (QC) were refrozen and stored at -80°C for at least 48 hours before being reanalyzed. When samples were reanalyzed due to QC failure, all within-person longitudinal samples were reanalyzed as well, on the second freeze/thaw cycle.

Dataset Information

This methods document applies to the following dataset(s) available from the ADNI repository:

| Dataset Name | Date Submitted |
|---|-----------------------|
| Fagan Lab – Longitudinal VILIP-1 CSF Analysis | 27 February 2015 |

References

1. Hulette CM, Welsh-Bohmer KA, Murray MG, et al. Neuropathological and neuropsychological changes in "normal" aging: evidence for preclinical Alzheimer disease in cognitively normal individuals. *J Neuropathol Exp Neurol*. Dec 1998;57(12):1168-1174.
2. Morris JC, Price JL. Pathologic correlates of nondemented aging, mild cognitive impairment, and early-stage Alzheimer's disease. *J Mol Neurosci*. Oct 2001;17(2):101-118.
3. Bateman RJ, Xiong C, Benzinger TL, et al. Clinical and biomarker changes in dominantly inherited Alzheimer's disease. *N Engl J Med*. Aug 30 2012;367(9):795-804.
4. Fagan AM, Xiong C, Jasielec MS, et al. Longitudinal change in CSF biomarkers in autosomal-dominant Alzheimer's disease. *Sci Transl Med*. Mar 5 2014;6(226):226ra230.
5. Villemagne VL, Burnham S, Bourgeat P, et al. Amyloid beta deposition, neurodegeneration, and cognitive decline in sporadic Alzheimer's disease: a prospective cohort study. *Lancet Neurol*. Apr 2013;12(4):357-367.
6. Tarawneh R, D'Angelo G, Macy E, et al. Visinin-like protein-1: diagnostic and prognostic biomarker in Alzheimer disease. *Ann Neurol*. Aug 2011;70(2):274-285.
7. Tarawneh R, Lee JM, Ladenson JH, Morris JC, Holtzman DM. CSF VILIP-1 predicts rates of cognitive decline in early Alzheimer disease. *Neurology*. Mar 6 2012;78(10):709-719.

About the Authors

This document was prepared by Dan Crimmins, Elizabeth Macy and Courtney Sutphen, Washington University in St. Louis, Departments of Pathology and Neurology. For more information please contact Courtney Sutphen at (314) 747-8396 or by email at sutphenc@neuro.wustl.edu.

Notice: This document is presented by the author(s) as a service to ADNI data users. However, users should be aware that no formal review process has vetted this document and that ADNI cannot guarantee the accuracy or utility of this document.